


LeviKam Active Sling

LAC-XS/S/M/L/XL-1



The LeviKam Active sling is compatible with hoists for people with disabilities. It features high-quality workmanship and a well-thought-out design. It is made of durable, safe polyester mesh that is highly water-permeable and dries quickly. The LeviKam Active sling allows the user to be transported in a seated position while providing back support. The user's position can be adjusted to their needs thanks to the four-point, adjustable suspension. A sewn-in maneuvering handle allows the user to be rotated while hanging.



LeviKam Active	XS size	S size	M size	L size	XL size
A size	84	86	96	107	122
B size	92	94	96	98	115
C size	22	24	26	28	31

*The values given are approximate due to differences in body structure.

Main advantages of vests:

- Maximum load capacity: 200 kg,
- Multi-level adjustment of the sling's hooks,
- Wide range of sizes available: from XS to XL,
- Slings feature a physiological cutout to facilitate toileting and hygiene procedures,
- Seats are available in two versions: a standard bath seat (made of polyester mesh) or a special comfort seat (made of solid polyester material, sewn upon request).

Stiffened head support L-UP-1

A practical and comfortable head stabilizing insert, designed to complement LeviKam vests, which do not normally have head support.



Instructional videos

Putting on slings in a sitting position:



Putting on slings in a lying position:

