



LeviKam Solid Sling

LSO-XS/S/M/L/XL-1



The LeviKam Solid sling is compatible with lifts for people with disabilities. It features high-quality workmanship and a well-thought-out design. It is made of durable, safe polyester mesh that is highly water-permeable and dries quickly. The LeviKam Solid sling allows the user to be transported in a seated position while providing back support. The user's position can be adjusted to their needs thanks to the four-point, adjustable suspension. A sewn-in maneuvering handle allows the user to be rotated while hanging.



LeviKam Solid	XS size	S size	M size	L size	XL size
A size	69	78	90	110	115
B size	81	88	96	97	104
C size	130	143	155	181	196
user weight	<20 kg	20-40 kg	40-70 kg	70-100 kg	>100 kg

*The values given are approximate due to differences in body structure.

Main advantages of vests:

- Maximum load capacity: 200 kg,
- Multi-level adjustment of the sling's hooks,
- Wide range of sizes available: from XS to XL,
- Slings feature a physiological cutout to facilitate toileting and hygiene procedures,
- Seats are available in two versions: a standard bath seat (made of polyester mesh) or a special comfort seat (made of solid polyester material, sewn upon request).

Instructional videos:

Putting on vests in a sitting position:



Putting on vests in a lying position:

